

## BT-1xxxSY/LN/Sci Brain slice chambers

The **BT-1xxxSY** brain slice chambers fit into the **Siskiyou Inc. 8090C** stage as well as our **MSAS-xxx** Stage adapters. Some are also available for the **Luigs & Neumann (LN)** stage and **Scientifica (Sci)** stages. These slice chambers all have inflow and outflow via a slit to promote laminar flow for flow rates up to 2ml/min (4ml/min for larger chambers). They come either with a smooth top or with a machined recess allowing electrode access at low angles.

Typically a transparent ITO heater (eg. **HI-25Dp**) will be used to directly heat the chamber from below. Perfusing solution is heated with the **HPRE2** Pre-heater connected to the inlet tube. Chambers have a thermistor groove to allow permanent placement of the thermistor sensor for measurement/control of the chamber temperature. Ground wires can either be placed in the outflow area or in the main working area to the side of the slice.

### Superior flow pattern

Turbulence in solution flowing through chambers can create regions with non-uniform drug concentrations. Our chamber designs demonstrate a predictable laminar flow pattern as in the example below.




Fig 1. shows dye containing solution flowing into **BT112B-45SY** fitted with a Warner Instruments, Inc 12mm harp (eg. SHD26H/10). Typical laminar flow pattern for solution front.


### Customization

Please call with any questions about modifying the chambers. We can for example vary the angle of the bevelled area to suit your lens or for example adapt a chamber for use with another microscope stage (see Fig. 2). With our advanced CAD/CAM software we have designed many chambers and can usually modify current chambers at no extra cost.



Fig 2. **BT-1-TB** chamber adapted to Luigs & Neumann stage.


## BT112B45/SY/Sci/LN



**Application:** brain slices-low electrode angles

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 12.7mm dia; 1mm lower level; 3.2mm height; 42° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**
- compatible with **12mm** harp (Warner Instruments, Inc SHD-26/xx)


## BT112B37HSY



**Application:** brain slices

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 12.7mm dia; 1mm lower level; 5mm height; 37° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**
- compatible with **12mm** harp (Warner Instruments, Inc SHD-26/xx)


## BT113B45/SY/Sci/LN



**Application:** brain slices-low electrode angles, 12mm cover slides

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 14.3mm diam; 1mm bottom height; 3.2mm height; 45° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**


## BT113B45HSY



**Application:** brain slices-12mm cover slides

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 14.3mm diam; 5mm height; 45° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**


## BT115B45/SY/Sci/LN



**Application:** brain slices-low electrode angles,

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 15.2mm diam; 3.2mm height; 45° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**
- compatible with **15mm** harp (Warner Instruments, Inc SHD-26G/xx)

## BT115B45HSY



**Application:** brain slices

- size platform 7.4 x 4.9 x 0.5cm
- size central chamber: 15.2mm diam; 5mm height; 45° bevel
- outflow chamber, ground & thermistor groove
- constructed of polycarbonate
- compatible with **HI-24p, HI-25Dp**
- compatible with **15mm** harp (Warner Instruments, Inc SHD-26G/xx)