


Cell MicroControls
P0 Box 11387, Norfolk, VA 23517
Tel: (800) 398-4262; (757) 622-0261
Fax: (757) 622-0262
WWW: <http://www.cellmc.com>; Email: info@cellmc.com
Equipment for cellular & electrophysiology research

Microscope stage adapters (partial list)

MSA-xxx Microscope stage adapters

The **MSA-xxx** microscope stage adapters allow hold the **BT-1-xxx** tissue chambers and **HLD-CS, HLD RG** culture chamber holders as well as the **BT-I55D** ITO stage. If you need an adapter for a different microscope please call us and we will make one for the same cost as a standard stage adapter.



MSA-MZ Marzhauser stage



MSA-NIKD Nikon Diaphot 10.8cm



MSA-OP Nikon Optiphot



MSA-OL11 Olympus 11cm



MSA-ZAX Zeiss Axiovert rect



MSA-TS100 Nikon




MSA-PR Prior automated stage



MSA-LEI/DM Leica DMIL



MSA-LEI/RD Leica round



MSA-ZFS Zeiss Axioskop FS



MSA-MECH1 Mechanical stage