



**Cell Micro Controls**  
 Norfolk, VA 23509  
 Tel: (757) 622-0261  
 Fax: (757) 622-0262  
 WWW: <http://www.cellmc.com>; Email: [info@cellmc.com](mailto:info@cellmc.com)  
*Equipment for cellular & electrophysiology research*

# Microscope stage adapters (partial list)

## MSA[S]-xxx Microscope stage adapters

The **MSA-xxx** microscope stage adapters hold the **BT-1-xxx** tissue chambers and **HLD-CS, HLD RG** culture chamber holders as well as the **BT-I55D** ITO stage. If you need an adapter for a different microscope please call us and we will make one for the same cost as a standard stage adapter. Adapters for chambers with the rectangular Siskiyou Inc outline (eg. **BT112B37HSY**) are designated **MSAS-xxx**.



**MSA-MZ Marzhauser stage**



**MSA-OP Nikon Optiphot**



**MSA-ZAX Zeiss Axiovert rect**



**MSA-PR Prior automated stage**



**MSA-LEI/RD Leica round**



**MSA-MECH1 Mechanical stage**



**MSA-NIKD Nikon Diaphot 10.8cm**



**MSA-OL11 Olympus 11cm**



**MSA-TS100 Nikon**



**MSA-LEI/DM Leica DMIL**



**MSA-ZFS Zeiss Axioskop FS**